

# EEGLAB- Quick Tutorial

Assaf Peleg, Dor Zazon

November 22, 2020

## First steps

First, download **EEGLAB** and acquire a relatively new version of **MATLAB**. You can download **EEGLAB** using the following link: <https://sccn.ucsd.edu/eeglab/download.php>. To access **EEGLAB** wiki go to: <https://sccn.ucsd.edu/wiki/EEGLAB>.

Next, open **MATLAB**, set your workspace to the folder in which you installed/extracted **EEGLAB** and write in the **Command Window** "eeglab" and press Enter.

## Examining a dataset

- 1) **Loading a dataset:** The first thing you might want to do is load a dataset, examine it's properties and take a look at it. To do so, click on **File** and then choose **Load existing dataset**. See **Figure 1**.




Figure 1: Loading an existing dataset

- 2) Next, the dataset will be loaded and its properties will be printed to the screen, e.g. the dataset that we have loaded has 32 channels(electrodes) and 154 events which are markers to align data with specific events that have occurred such as stimulus on the screen. See **Figure 2**.




Figure 2: Dataset properties

- 3) **Looking at the data:** The following thing you might want to do is to look at the data with your own eyes. To do so, click on **Plot** and then choose **Channel data (scroll)**. See **Figure 3**.




Figure 3: Looking at the data

- 4) A window will open. On the left you'll see the names of the channels/electrodes, in the center you'll see the actual data and on the bottom you'll see arrows allowing you to scroll it. You'll also notice colored lines- these are the events we mentioned earlier. See **Figure 4**.




Figure 4: Scrolling the data

## Cleaning the data

- 1) **Removing bad sections:** A look at the data might reveal some bad and noisy sections(in all the electrodes at the same time), and we might want to manually delete them(This will affect the phase of the waves because we cut the data but it will not be significant). This can be done by selecting them using the mouse(click on the start of the section, hold your mouse and drag it to the end of the section). Then, click on the **Reject** button. See **Figure 5**.




Figure 5: Rejecting bad sections

- 2) Next the scrolling window will close and a new window will open. This window will allow you to choose if to overwrite your current dataset with the changes by selecting overwrite and ok or create a new one by clicking ok without selecting anything(See **Figure 6**)- **Explanation:** When working in **EEGLAB**, datasets can be loaded and used in parallel. Each time you preform an action, it will be preformed on the current dataset. To choose a dataset go to **Datasets**. For example, we chose not to overwrite our datasets so now we can choose between two different datasets. See **Figure 7**.




Figure 6: Saving the changes




Figure 7: Choosing a dataset

- 3) **Removing bad electrodes:** Our next order of business is to look for bad electrodes as some of them might be noisy and full of garbage altogether. After we scroll the data once more, and decide that there is such an electrode, we would want to delete it. To do so, click on **Edit** and then choose **Select data**. See **Figure 8**.




Figure 8: Selecting data to reject

- 4) Next, check the box of **"on→remove these"** for **Channel range** and click on the button with the ellipses(...) to choose the electrodes you want to delete. Finally click on ok and choose if to overwrite the dataset or create a new one. See **Figure 9** and **Figure 10**.




Figure 9: After checking the box and clicking on the ellipses




Figure 10: After picking an electrode

- 5) If you scroll the data now, you won't see that electrode anymore. See **Figure 11**.




Figure 11: F3 electrode is now gone

## Filtering the data

- 1) **HIGH PASS Filter**: Next, we'll want to filter the data by searching for which frequencies we have undesirable peaks and checking if they need to be filtered or not (We'll search for the frequencies which we would want to always remove). However, first thing first and we'll need to look at the data again. This time will use a PSD graph. To do so, click on **Plot** and then choose **Channel spectra and maps**. See **Figure 12**.




Figure 12: Channel spectra and maps

- 2) A window will open which will allow us to choose several options for the PSD graph. See **Figure 13**.




Figure 13: PSD graph options

- 3) For example, we would like to view the PSD graph for frequencies between 0hz and 64hz, and therefore we change the plotting range and click ok. See **Figure 14**.




Figure 14: PSD graph plotting range- 0hz to 64hz

- 4) Next, the PSD graph will open. First, we can see the alpha waves at 10hz and we see an interesting anomaly at 60hz. This is the American power supply frequency where in Israel it will be 50hz. See **Figure 15**. We would surely want to filter that. We can do this using a **HIGH PASS Filter**.




Figure 15: PSD



- 5) To do so, go to **Tools**, **Filter the data** and then **Basic FIR**. See **Figure 16**.




Figure 16: Basic FIR

- 6) The filter we'll want to do mostly is 0.5hz to 40hz as we are not interested in higher or lower frequencies. To do so, enter 0.5 to the low filter and 40 to the high filter and click ok. Click ok to create a new dataset as always when the saving window opens. See **Figure 17**.




Figure 17: 0.5hz to 40hz HIGH PASS filter

- 7) Next, a window will open plotting the filter being used. See **Figure 18**. However, to see the actual affect on our data we'll need to plot another PSD graph (Don't forget to change the plotting range as before). See **Figure 19**. As can be seen from the PSD graph, the 60hz peak was removed.




Figure 18: Filter plot




Figure 19: PSD graph

- 8) **Notch Filter**: Another way to filter is to use a notch filter. A notch filter will allow us to discard a range of frequencies. First, we'll return to the previous dataset before the affects of the high pass filter. In order to do so, go to **Datasets** and then choose dataset 3. See **Figure 20**.




Figure 20: Choosing the previous dataset

- 9) Next, to apply the filter go to **Tools**, **Filter the data** and then **Basic FIR**. See **Figure 21**.




Figure 21: Opening the filter window again

- 10) Then, check the box near the notch filter and enter in the low filter input box 57 and in the high filter input box enter 63. See **Figure 22**.




Figure 22: Applying a notch filter

- 11) Next, a window will open plotting the filter being used. However, to see the actual affect on our data we'll need to plot another PSD graph(Don't forget to change the plotting range as before). See **Figure 23**. As can be seen from the PSD graph, the 60hz peak was removed, but this time the filter was narrower and replaced the 60hz peak with a negative one.




Figure 23: PSD graph

- 12) Before we continue, let's return to dataset 4 in which we preformed the high pass filter.

## Final touches

- 1) **Average reference:** The following thing we might want to do is removing relatively robustly bias and noise that affects all the electrodes altogether such as eye blinks. To do so, we can calculate the average reference and remove it. Click on **Tools** and choose **Re-reference the data**. See **Figure 24**.




Figure 24: Re-referencing the data

- 2) Make sure **Compute the average reference** is selected, click ok and save the new dataset. See **Figure 25**.




Figure 25: Average reference

- 3) **Interpolation:** Next, after we cleaned the data with various methods, we need to handle our deleted electrode. Recall that we removed F3 electrode as we thought it was too noisy and full of garbage data. However, if we have many test subjects, and in each subject we delete a different electrode, that might cause problems. Therefore we cannot just remove noisy electrodes. To mitigate this we can use interpolation. First, to make this clear, choose the first dataset and then click on **Plot, Channel locations** and then choose **By name**. See **Figure 26**.




Figure 26: Plotting channel locations

- 4) The way we will mitigate our issue is to use nearby electrodes data from FC5, FC1, Fz, etc... to figure out what might the F3 electrode would have reported. See **Figure 27**.




Figure 27: Channel location plot

- 5) So, we figured out we need interpolation. First, choose our latest dataset number 6. Now, click on **Tools** and then choose **Interpolate electrodes**. See **Figure 28**.




Figure 28: Interpolate electrodes

- 6) Next, the interpolation needs to know the location of the F3 electrode. However, as we deleted it, we don't know its location. The solution to this is to use our other datasets. For example, our first dataset must contain F3 location. So, in the window that opens choose **Use all channels from other dataset**, in the window that opens enter 1 and click ok. See **Figure 29**.




Figure 29: Use all channels from other dataset

- 7) If all went well you should see at the top of the window that the F3 electrode was figured out automatically as the one to be interpolated. See **Figure 30**.




Figure 30: Interpolation

- 8) Now click ok and choose not to overwrite the dataset. If you plot now the data you should see that the F3 electrode is back. See **Figure 31**.




Figure 31: F3 electrode is back

- 9) **Re-sampling:** A final touch might be to re-sample to data in order to save space and later in time complexity. You could do this by going to **Tools** and choosing **Change sampling rate**.



## ICA, RANSAC and ASR

- 1) **ICA**: After we cleaned our data, we can use several algorithms to utilize it(ICA) or continue cleaning it(ASR, RANSAC). One such algorithm is ICA. ICA allows us separate our data from all the electrodes to different independent components such as eye blinks. For more info: [https://scn.ucsd.edu/wiki/Chapter\\_09:\\_Decomposing\\_Data\\_Using\\_ICA](https://scn.ucsd.edu/wiki/Chapter_09:_Decomposing_Data_Using_ICA)
- 2) To perform ICA, go to **Tools** and then choose **Decompose data by ICA**. See **Figure 32**.




Figure 32: Decompose data by ICA

- 3) In the window that opens click ok(this might take some time so be patient, you can see the progress of the ICA training in the **MATLAB Command Window**). See **Figure 33**.




Figure 33: ICA

- 4) Next, we can plot the different ICA components using a scroll or components maps. See **Figure 34** and **Figure 35**.




Figure 34: Scroll




Figure 35: Components maps

- 5) In the scroll plot, we can see at the left that now instead of electrodes we have the ICA components by number. See **Figure 36**.




Figure 36: Scroll plot

- 6) In the components maps plot we can see a map plot of each components. See **Figure 37**. For example, we can see the component 20 is hot in the front of the face, this might be from eye blinks or from other muscles of the face- therefore we might want to delete that component. For more info see the wiki.




Figure 37: Components maps plot

- 7) **ASR & RANSAC**: The ASR and RANSAC algorithms might help you automatically clean noisy data. For more info, see the wiki.

## Getting the code & Saving the dataset

- 1) **Exporting to code:** To export your actions into code, you can write "eegh" in the **Command Window** of **MATLAB**. You might want to preform the actions without plotting as these actions will be exported as well.
- 2) **Saving the dataset:** After you are finished, you can save your dataset. To save your dataset, go to **File** and choose one of several save options. See **Figure 38**.




Figure 38: Saving the dataset

Good luck and enjoy from Assaf & Dor